

2016-07-14

EST Medium

Pers. comm. Mix
Modified by MZCH

Composition of EST Medium.

| | stock solutions (1 L) | | volume of stock for 1 L nutrient solution | final concentration |
|--|--------------------------|--------|--|------------------------|
| | [mmol/L] | [g/L] | [mL] | [mmol/L] |
| KNO ₃ | 1000 | 101.11 | 2 | 2 |
| CaCl ₂ ·6H ₂ O | 22.82 | 5 | 0.1 | 0.002 |
| FeCl ₃ | 6.16 | 1 | 1.5 | 0.009 |
| MgSO ₄ ·7H ₂ O | 40.57 | 10 | 1 | 0.04 |
| (NH ₄) ₂ HPO ₄ | 151.45 | 20 | 2 | 0.303 |
| Na ₂ EDTA·2H ₂ O (Titrplex® III) | 26.86 | 10 | 1 | 0.027 |
| micronutrient solution ¹ | | | 2.5 | |
| soil extract ² | | | 100 | |
| NP-H ₂ O ³ | | | 889.9 | |

Adjust pH to 6.

¹micronutrient solution stock (400x, 1 L)

| | stock [mmol/L] | stock [g/L] |
|--|----------------|-------------|
| CuSO ₄ ·5H ₂ O | 0.016 | 0.004 |
| Co(NO ₃) ₂ ·6H ₂ O | 0.137 | 0.04 |
| H ₃ BO ₃ | 6.47 | 0.4 |
| MnSO ₄ ·4H ₂ O | 0.36 | 0.08 |
| Na ₂ MoO ₄ ·2H ₂ O | 0.165 | 0.04 |
| ZnSO ₄ ·7H ₂ O | 1.39 | 0.4 |
| ddH ₂ O ⁴ , ad 1 L | | |

²preparation of soil extract:

Weigh 50 g of beech forest soil in a 1 L Erlenmeyer flask. Add 625 mL of dest. H₂O. Heat the soil solution and keep at 100 °C for 5 minutes. Place a piece of pleated filter paper on top of a second 1 L Erlenmeyer flask and fill a spatula tip of CaCO₃ on the filter. Filter the soil extract (if necessary overnight) and then stir for 15 min. The extract is then centrifuged for 15 minutes at 2500 g and 20 °C. Transfer and aliquot the supernatant into 50 mL Falcon tubes. Store at -20°C.

The volume of the added soil extract should be adapted to the nutrient requirements of the cultured strains.

³NP-H₂O nanopure water, Purelab Pulse (ELGA Lab water, Celle, Germany)

⁴ddH₂O double distilled water

Reference

Pers. Comm. Marianne Mix (Univ. Hamburg).